



model 7554

Dust-Proof Personal Emergency Eyewash Station

FEATURES & BENEFITS

DUST PROTECTION

Tightly-sealed durable enclosure protects eyewashes from dust and debris, keeping bottles clean and hygienic in dirty environments.

IMMEDIATE FLUSHING

The eyewashes provide immediate flushing fluid while the injured person is en route to a primary emergency station.

DUAL EYE PROTECTION

With two ergonomic streams, the DUO eyewash provides simultaneous eye coverage allowing both eyes to be treated with an advanced rinse process for supplementary flushing support.

STERILE SALINE EFFECTIVENESS

The sterile sodium chloride Saline solution rinses the eyes to reduce the risk of severe injury caused by chemicals or particulate contamination.

PH NEUTRAL ADVANTAGE


A sterile pH Neutral phosphate solution is able to neutralize any acidic or alkali substance that enters the eyes to help restore a proper pH balance (only available pH neutralizer on the market). Recommended for chemicals with Ph levels above 10 or below 4, though it can also be used for pH neutral contaminants.

AS4775 and ANSI COMPLIANCE

AS4775 and ANSI compliant for supplemental equipment.

MOUNTING

Wall-mounted design allows for flexible placement next to the potential hazard. Includes two wall hooks for easy installation, and wall-mounted sign with universal language for placement in any situation.



SPECIFICATIONS

Model 7554 is a wall-mounted, dust-proof, dual bottle personal eyewash station that includes one 500 ml (16.9 fl oz) Sterile pH Neutral solution bottle and one 1000 ml (33.8 fl oz) Sterile Saline solution bottle in a tightly-sealed enclosure. Recommended wash time: 2 minutes with the pH Neutral solution and 5 minutes with the Saline solution. Includes a visual assistance mirror. Fluid carries a 3 year shelf life. Store at room temperature 5-35° C (41-95° F)

APPLICATIONS

Haws Personal Emergency Eyewashes are for use when dirt, dust or other contaminants come into contact with the eyes. Ideal for quick access to first aid in case of eye injury. Mounted at the immediate point of need, this personal eyewash will provide immediate flushing until a primary station can be reached.

OPTIONS

- ❑ Personal Emergency Eyewash (Replacement Bottle): Model 7583 is a replacement bottle of the 500 ml (16.9 fl oz) sterile pH Neutral solution. With two ergonomic streams, the DUO eyewash provides simultaneous eye coverage allowing for both eyes to be treated with an advanced rinse process for supplementary flushing support. The sterile pH Neutral phosphate solution is able to neutralize any acidic or alkali substance that enters the eyes to help restore a proper pH balance.
- ❑ Personal Emergency Eyewash (Replacement Bottle): Model 7584 is a replacement bottle of the 1000 ml (33.8 fl oz) sterile Saline solution. With two ergonomic streams, the DUO eyewash provides simultaneous eye coverage allowing for both eyes to be treated with an advanced rinse process for supplementary flushing support.

For more information, visit www.allvalve.com.au or call +61 2 8543 9811